

# Mouse Lung Digestion to Single-Cell Suspension

Abdullah Lab, IMMEI, University Hospital Bonn

2025-11-20

**Protocol ID:** MUS-LNG-DIG-001

**Version:** v1.0

**Author:** Dillon Corvino

## Purpose

This protocol describes the enzymatic and mechanical dissociation of mouse lung tissue to obtain single-cell suspensions suitable for flow cytometry, FACS sorting, and single-cell sequencing. Lung tissue is fibrous and elastic; therefore, digestion requires both **collagenase** and **DNase** together with thorough mechanical disruption.

## Critical notes (read before starting)

- Keep all buffers and tubes **cold** except during digestion at **37 °C**.
- Use **fresh, pre-warmed digestion buffer** with DNase added on the day of use.
- Lung tissue contains **trapped blood**; incomplete perfusion increases RBC contamination.
- Do not over-digest: limit digestion to **30–35 minutes**, as prolonged incubation reduces viability.
- After digestion, perform all downstream steps on **ice**.
- Use **70 µm strainers** and do not reuse between lungs.

## Approximate timing

- Dissection: **5 min**
- Mince + pre-digestion handling: **5 min**
- Enzymatic digestion: **30–35 min**
- Mechanical processing and washes: **20–25 min**

Total hands-on time: **~60 minutes per lung** (less when batching).

## Table of contents

<b>Purpose</b>	<b>1</b>
<b>Critical notes (read before starting)</b>	<b>1</b>
<b>Approximate timing</b>	<b>1</b>
<b>Procedure</b>	<b>3</b>
STEP 1 – Lung collection . . . . .	3
STEP 2 – Tissue mincing . . . . .	3
STEP 3 – Enzymatic digestion . . . . .	3
STEP 4 – Mechanical dissociation . . . . .	3
STEP 5 – RBC removal (if required) . . . . .	4
STEP 6 – Final wash and resuspension . . . . .	4
<b>Buffers used</b>	<b>5</b>
<b>Materials</b>	<b>6</b>
Reagents . . . . .	6
Disposables . . . . .	6
Equipment . . . . .	6
<b>Troubleshooting</b>	<b>7</b>
<b>Safety (brief)</b>	<b>8</b>
<b>Version history</b>	<b>8</b>

## Procedure

### STEP 1 – Lung collection

1. Euthanise the mouse using an approved method.
2. Open the thoracic cavity to expose the lungs.
3. (Optional but recommended) Perfuse via the right ventricle with **cold PBS** until lungs visibly blanch, reducing RBC content.
4. Excise the entire lung block and place it in a **10 cm petri dish** containing **cold PBS or FACS/MACS Buffer**.

#### Tip

Perfusion markedly improves final leukocyte purity and reduces RBC lysis requirements.

### STEP 2 – Tissue mincing

1. Blot the lungs briefly on sterile tissue to remove excess fluid.
2. Transfer to a fresh **10 cm dish** and mince thoroughly using scissors or scalpels until small fragments are obtained.
3. Transfer the minced lung into a **50 mL Falcon tube** containing **5 mL Lung Digestion Buffer** (see buffer protocol BUF-LNG-DIG-v1.0).

### STEP 3 – Enzymatic digestion

1. Place tubes on a **37 °C benchtop shaker** at **250 rpm**.
2. Incubate for **30–35 minutes**.

#### Warning

Over-digestion can decrease viability and negatively impact epithelial and myeloid populations.

### STEP 4 – Mechanical dissociation

1. After digestion, add **FACS/MACS Buffer** to bring total volume to **25 mL**.
2. Pour the suspension onto a **70 µm cell strainer** placed on a **10 cm dish**.
3. Mechanically mash the digested lung through the mesh using the **back of a 5 mL syringe plunger**.
4. Rinse the strainer with **cold FACS/MACS Buffer** until no visible tissue remains.
5. Transfer the filtered cell suspension back into a **50 mL tube**.
6. Top up to **50 mL** with cold FACS/MACS Buffer.
7. Centrifuge at **400 g for 5 minutes at 4 °C**.
8. Carefully discard the supernatant.

## STEP 5 – RBC removal (if required)

Lung digests vary in RBC content. If erythrocytes are visible:

1. Add **1 mL 1× RBC Lysis Buffer** (e.g. HybriMax) per lung.
2. Mix gently and incubate for **3–5 minutes at room temperature**.
3. Immediately top up to **50 mL** with cold FACS/MACS Buffer.
4. Centrifuge at **400 g for 5 minutes at 4 °C**.
5. Discard supernatant.

### **i** Note

If the lung was perfused effectively, this step may not be needed.

## STEP 6 – Final wash and resuspension

1. Loosen pellet by gentle flicking or brief vortexing.
2. Resuspend in **5 mL FACS/MACS Buffer** (or other desired volume).
3. Filter again through a **fresh 70 µm strainer** if clumping is observed.
4. Keep the suspension **on ice** until staining, sorting, or single-cell workflows.

## Buffers used

- **Lung Digestion Buffer** – BUF-LNG-DIG-v1.0
  - Typically DMEM (or RPMI) + Collagenase D or Collagenase IV + DNase I.
- **FACS/MACS Buffer** – BUF-FACS-v1.0
- **RBC Lysis Buffer** – BUF-RBC-HYB-v1.0

## Materials

### Reagents

Reagent	Supplier / Cat#	Notes
Lung Digestion Buffer	In-house	Pre-warmed to 37 °C
FACS/MACS Buffer	In-house	Cold
RBC Lysis Buffer (1×)	[TBD]	Use only if needed
PBS (1×)	Any	For perfusion and rinsing

### Disposables

Item	Notes
50 mL Falcon tubes	1 per lung
10 cm petri dishes	For mincing & filtering
70 µm strainers	Do not reuse
5 mL syringes (no needle)	For mechanical dissociation
Scissors, scalpels	For mincing
Transfer pipettes	For washing & transfers

### Equipment

Equipment	Notes
Benchtop shaker/incubator	37 °C, 250 rpm
Benchtop centrifuge	400 g at 4 °C
Biosafety cabinet (Class II)	Recommended
Timer	For digestion and lysis steps

## Troubleshooting

Issue	Possible cause	Suggested solution
Persistent clumps	Incomplete digestion	Increase mechanical dissociation; do not over-digest
High RBC contamination	Lungs not perfused	Perform RBC lysis or perfuse before collection
Low viability	Over-digestion or RT handling	Limit digestion to 30–35 min; keep washes cold
Excess debris	Harsh mechanical disruption	Use gentle, repetitive plunger strokes

## Safety (brief)

- Treat mouse tissues as potentially infectious and process according to institutional S1/S2 rules.
- Enzyme mixtures (Collagenase, DNase) may cause irritation; wear gloves and eye protection.
- Dispose of sharps, tissues, and enzymatic waste according to facility guidelines.

## Version history

Version	Date	Author	Change summary
v1.0	2025-11-20	Dillon Corvino	Initial upgraded version for lung digestion