

# K-562 Cell Line Culture Protocol

Abdullah Lab, IMMEI, University Hospital Bonn

**Protocol ID:** HUM-K562-CULTURE-001

**Version:** v1.0

**Author:** Dillon Corvino

**Date:** 2025-11-21

## Purpose

This protocol describes the thawing, routine culture, maintenance, passaging, and cryopreservation of the human K-562 erythroleukemia cell line. K-562 cells grow in suspension and require regular monitoring of density to maintain viability and consistent growth rates.

## Critical notes (read before starting)

- K-562 cells grow as **suspension cultures** and should remain in exponential growth phase.
- Do not allow cultures to exceed recommended density, as nutrient depletion reduces viability and alters phenotype.
- K-562 cells are of human origin—treat all cultures and waste as BSL2.
- Always pre-warm media and reagents to **37 °C** before use.
- Ensure gentle handling; harsh centrifugation or over-pipetting can damage suspension cells.
- Cryopreservation requires **10% DMSO**, and cells should be frozen using a controlled-rate freezing container.

### Warning

Overgrowth (tight clumps, very dense cultures) leads to poor viability and altered functional responses. Split cultures **before** they exceed recommended density.

### Tip

Record passage number and split ratios at every passage. This improves reproducibility when using K-562 cells as targets or stimulator cells.

## **Approximate timing**

- Thawing: 15–20 min
- Routine medium change / split: 10–15 min
- Cryopreservation: 20–30 min
- Total hands-on time per session: ~15–30 min

## Table of contents

<b>Purpose</b>	<b>1</b>
<b>Critical notes (read before starting)</b>	<b>1</b>
<b>Approximate timing</b>	<b>2</b>
<b>Procedure</b>	<b>4</b>
Step 1 – Reagents and medium preparation . . . . .	4
Step 2 – Thawing K-562 cells . . . . .	4
Step 3 – Routine culture and passaging . . . . .	4
Recommended density . . . . .	4
Procedure . . . . .	5
Step 4 – Short-term storage / holding conditions . . . . .	5
Step 5 – Cryopreservation of K-562 cells . . . . .	5
Preparation of freezing medium . . . . .	5
Freezing procedure . . . . .	5
<b>Materials</b>	<b>6</b>
Reagents . . . . .	6
Disposables . . . . .	6
Equipment . . . . .	6
<b>Buffers used</b>	<b>6</b>
<b>Troubleshooting</b>	<b>6</b>
<b>Safety (brief)</b>	<b>7</b>
<b>Version history</b>	<b>7</b>

## Procedure

### Step 1 – Reagents and medium preparation

1. Pre-warm complete culture medium:
  - RPMI 1640
  - 10% FCS
  - 1% Pen/Strep
2. Pre-warm PBS (if washing is needed).
3. Label flasks with:
  - Cell line name (K-562)
  - Passage number
  - Date
  - Operator initials

### Step 2 – Thawing K-562 cells

1. Pre-warm **10–15 mL** of complete medium to 37 °C.
2. Thaw a cryovial of K-562 cells in a 37 °C water bath:
  - Gently swirl and remove when only a small ice crystal remains.
3. Transfer contents into a 15 mL tube containing **10 mL** pre-warmed complete medium.
4. Centrifuge at **300 g for 5 min** at room temperature.
5. Remove supernatant and gently resuspend the pellet in **10–15 mL** fresh complete medium.
6. Transfer cells to a T25 or T75 flask and incubate at **37 °C, 5% CO<sub>2</sub>**.

#### **i** Note

The centrifugation step removes DMSO, preventing toxicity during recovery.

### Step 3 – Routine culture and passaging

#### Recommended density

- Maintain between  **$2 \times 10^5$**  and  **$1 \times 10^6$**  cells/mL.

## Procedure

1. Gently resuspend the culture by pipetting up and down.
2. Count cells using trypan blue or an automated counter.
3. If density  $> 1 \times 10^6$  cells/mL:
  - Split culture **1:2 to 1:5** depending on growth.
4. Add fresh, pre-warmed complete medium to adjust culture to:
  - **$2\text{--}3 \times 10^6$  cells/mL**
5. Return to incubator (37 °C, 5% CO<sub>2</sub>).

### Warning

Do not allow cultures to sit for more than 3–4 days without splitting. Nutrient depletion leads to apoptosis and growth arrest.

## Step 4 – Short-term storage / holding conditions

- Cultures may be held at  **$2\text{--}4 \times 10^6$  cells/mL** for 24–48 h without splitting.
- If cultures appear slow or viability drops:
  - Increase medium volume
  - Spin down and resuspend in entirely fresh medium

## Step 5 – Cryopreservation of K-562 cells

### Preparation of freezing medium

- **90% FCS + 10% DMSO**  
(or use BUF-GEN-FREEZE-MIX-001 diluted 1:1 with cell suspension)

### Freezing procedure

1. Collect cells in exponential growth phase.
2. Centrifuge at **300 g for 5 min.**
3. Resuspend in freezing medium at  **$1\text{--}10 \times 10^6$  cells/mL.**
4. Aliquot **1 mL** into labelled cryovials.
5. Place vials into a **controlled-rate freezing container** (e.g. Mr Frosty).
6. Freeze at  **$-70/-80$  °C overnight.**
7. Transfer vials to **liquid nitrogen** for long-term storage.

## Materials

### Reagents

Reagent	Supplier	Cat. #	Notes
RPMI 1640	Various	—	Base medium
Fetal calf serum (FCS)	Various	—	10% v/v in medium
Penicillin/Streptomycin	Various	—	1% v/v final
DMSO (cell culture grade)	Sigma/others	—	10% final in freezing medium
PBS, 1×	Various	—	For washing (optional)

### Disposables

Item	Specification
T25 / T75 flasks	Vent-cap, sterile
15 mL / 50 mL tubes	Sterile
Cryovials	LN2-compatible
Pipette tips	Sterile (10 $\mu$ L – 1000 $\mu$ L)
Serological pipettes	Sterile (5–25 mL)

### Equipment

Equipment	Specification / Notes
CO incubator	37 °C, 5% CO
Benchtop centrifuge	300 g for washing
Biological safety cabinet	BSL2 handling
Water bath	37 °C for thawing
Cell counter/hemocytometer	For counting
Controlled-rate freezer	e.g. Mr Frosty
Liquid nitrogen storage	Long-term cryostorage

### Buffers used

- **BUF-GEN-FREEZE-MIX-001** – optional alternative for freezing
- PBS (standard laboratory reagent)

### Troubleshooting

Issue	Possible cause	Solution
Low viability after thaw	Slow thawing or DMSO exposure	Thaw rapidly; wash immediately after thaw
Clumping	High density or depleted medium	Dilute, split earlier, or replace medium entirely
Slow growth	Old medium, poor serum batch	Prepare fresh medium; test new batch of FCS
Sudden drop in viability	Contamination or overgrowth	Check culture; discard and restart from stock

## Safety (brief)

- Handle all human-derived cell lines using BSL2 precautions.
- Wear gloves, lab coat, and eye protection.
- Dispose of liquid and solid waste as biohazardous material.
- DMSO increases skin permeability—avoid contact.

## Version history

Version	Date	Author	Changes
v1.0	2025-11-21	Dillon Corvino	Initial Quarto protocol for K-562 culture