

K-562 Cell Line Culture Protocol

Abdullah Lab, IMMEI, University Hospital Bonn

Protocol ID: HUM-K562-CULTURE-001

Version: v1.0

Author: Dillon Corvino

Date: 2025-11-21

Purpose

This protocol describes the thawing, routine culture, maintenance, passaging, and cryopreservation of the human K-562 erythroleukemia cell line. K-562 cells grow in suspension and require regular monitoring of density to maintain viability and consistent growth rates.

Critical notes (read before starting)

- K-562 cells grow as **suspension cultures** and should remain in exponential growth phase.
- Do not allow cultures to exceed recommended density, as nutrient depletion reduces viability and alters phenotype.
- K-562 cells are of human origin—treat all cultures and waste as BSL2.
- Always pre-warm media and reagents to **37 °C** before use.
- Ensure gentle handling; harsh centrifugation or over-pipetting can damage suspension cells.
- Cryopreservation requires **10% DMSO**, and cells should be frozen using a controlled-rate freezing container.

Warning

Overgrowth (tight clumps, very dense cultures) leads to poor viability and altered functional responses. Split cultures **before** they exceed recommended density.

Tip

Record passage number and split ratios at every passage. This improves reproducibility when using K-562 cells as targets or stimulator cells.

Approximate timing

- Thawing: 15–20 min
- Routine medium change / split: 10–15 min
- Cryopreservation: 20–30 min
- Total hands-on time per session: ~15–30 min

Table of contents

Purpose	1
Critical notes (read before starting)	1
Approximate timing	2
Procedure	4
Step 1 – Reagents and medium preparation	4
Step 2 – Thawing K-562 cells	4
Step 3 – Routine culture and passaging	4
Recommended density	4
Procedure	5
Step 4 – Short-term storage / holding conditions	5
Step 5 – Cryopreservation of K-562 cells	5
Preparation of freezing medium	5
Freezing procedure	5
Materials	6
Reagents	6
Disposables	6
Equipment	6
Buffers used	6
Troubleshooting	6
Safety (brief)	7
Version history	7

Procedure

Step 1 – Reagents and medium preparation

1. Pre-warm complete culture medium:
 - RPMI 1640
 - 10% FCS
 - 1% Pen/Strep
2. Pre-warm PBS (if washing is needed).
3. Label flasks with:
 - Cell line name (K-562)
 - Passage number
 - Date
 - Operator initials

Step 2 – Thawing K-562 cells

1. Pre-warm **10–15 mL** of complete medium to 37 °C.
2. Thaw a cryovial of K-562 cells in a 37 °C water bath:
 - Gently swirl and remove when only a small ice crystal remains.
3. Transfer contents into a 15 mL tube containing **10 mL** pre-warmed complete medium.
4. Centrifuge at **300 g for 5 min** at room temperature.
5. Remove supernatant and gently resuspend the pellet in **10–15 mL** fresh complete medium.
6. Transfer cells to a T25 or T75 flask and incubate at **37 °C, 5% CO₂**.

i Note

The centrifugation step removes DMSO, preventing toxicity during recovery.

Step 3 – Routine culture and passaging

Recommended density

- Maintain between **2×10^5** and **1×10^6** cells/mL.

Procedure

1. Gently resuspend the culture by pipetting up and down.
2. Count cells using trypan blue or an automated counter.
3. If density $> 1 \times 10^6$ cells/mL:
 - Split culture **1:2 to 1:5** depending on growth.
4. Add fresh, pre-warmed complete medium to adjust culture to:
 - **$2\text{--}3 \times 10^6$ cells/mL**
5. Return to incubator (37 °C, 5% CO₂).

Warning

Do not allow cultures to sit for more than 3–4 days without splitting. Nutrient depletion leads to apoptosis and growth arrest.

Step 4 – Short-term storage / holding conditions

- Cultures may be held at **$2\text{--}4 \times 10^6$ cells/mL** for 24–48 h without splitting.
- If cultures appear slow or viability drops:
 - Increase medium volume
 - Spin down and resuspend in entirely fresh medium

Step 5 – Cryopreservation of K-562 cells

Preparation of freezing medium

- **90% FCS + 10% DMSO**
(or use BUF-GEN-FREEZE-MIX-001 diluted 1:1 with cell suspension)

Freezing procedure

1. Collect cells in exponential growth phase.
2. Centrifuge at **300 g for 5 min.**
3. Resuspend in freezing medium at **$1\text{--}10 \times 10^6$ cells/mL.**
4. Aliquot **1 mL** into labelled cryovials.
5. Place vials into a **controlled-rate freezing container** (e.g. Mr Frosty).
6. Freeze at **$-70/-80$ °C overnight.**
7. Transfer vials to **liquid nitrogen** for long-term storage.

Materials

Reagents

Reagent	Supplier	Cat. #	Notes
RPMI 1640	Various	—	Base medium
Fetal calf serum (FCS)	Various	—	10% v/v in medium
Penicillin/Streptomycin	Various	—	1% v/v final
DMSO (cell culture grade)	Sigma/others	—	10% final in freezing medium
PBS, 1×	Various	—	For washing (optional)

Disposables

Item	Specification
T25 / T75 flasks	Vent-cap, sterile
15 mL / 50 mL tubes	Sterile
Cryovials	LN2-compatible
Pipette tips	Sterile (10 μ L – 1000 μ L)
Serological pipettes	Sterile (5–25 mL)

Equipment

Equipment	Specification / Notes
CO incubator	37 °C, 5% CO
Benchtop centrifuge	300 g for washing
Biological safety cabinet	BSL2 handling
Water bath	37 °C for thawing
Cell counter/hemocytometer	For counting
Controlled-rate freezer	e.g. Mr Frosty
Liquid nitrogen storage	Long-term cryostorage

Buffers used

- **BUF-GEN-FREEZE-MIX-001** – optional alternative for freezing
- PBS (standard laboratory reagent)

Troubleshooting

Issue	Possible cause	Solution
Low viability after thaw	Slow thawing or DMSO exposure	Thaw rapidly; wash immediately after thaw
Clumping	High density or depleted medium	Dilute, split earlier, or replace medium entirely
Slow growth	Old medium, poor serum batch	Prepare fresh medium; test new batch of FCS
Sudden drop in viability	Contamination or overgrowth	Check culture; discard and restart from stock

Safety (brief)

- Handle all human-derived cell lines using BSL2 precautions.
- Wear gloves, lab coat, and eye protection.
- Dispose of liquid and solid waste as biohazardous material.
- DMSO increases skin permeability—avoid contact.

Version history

Version	Date	Author	Changes
v1.0	2025-11-21	Dillon Corvino	Initial Quarto protocol for K-562 culture