

# ACK Red Blood Cell Lysis Buffer

Abdullah Lab, IMMEI, University Hospital Bonn

2026-01-09

**Protocol ID:** BUF-ACK-LYSIS-001

**Version:** v1.1

**Author:** Dillon Corvino

## Purpose

ACK lysis buffer is used to selectively lyse erythrocytes in single-cell suspensions derived from mouse or human tissues (e.g. liver digests, spleen, PBMC preparations). Under brief incubation, the buffer efficiently removes red blood cells while preserving leukocyte viability and integrity.

## Working buffer composition

Component	Final concentration	Notes
Ammonium chloride (NH Cl)	X mM	Primary erythrocyte-lysing component
Potassium bicarbonate (KHCO )	X mM	Buffering capacity
EDTA	X mM	Chelates divalent cations; reduces cell clumping
Sterile distilled water	—	Solvent
pH	7.2	Adjusted after dissolution

### Note

This formulation reflects the standard ACK buffer used in immunology workflows, with explicit pH adjustment to ensure reproducibility across batches.

## Preparation

### Reagents

Reagent	Form	Notes
NH Cl	Powder	Analytical or cell-culture grade
KHCO	Powder	Hygroscopic; store tightly sealed
EDTA	Powder	Cell-culture grade
Distilled water	—	Sterile (bideist)
HCl	—	For pH adjustment

### Preparation of 1 L ACK lysis buffer

1. Use a clean **1 L glass bottle** or beaker.
2. Add the following dry reagents:
  - **8.26 g NH Cl**
  - **1.00 g KHCO**
  - **0.037 g (37 mg) EDTA**
3. Add **1 L sterile distilled (bideist) water**.
4. Mix thoroughly by stirring or gentle swirling until all components are fully dissolved.
5. Measure pH and **adjust to pH 7.2** using dilute HCl.
6. **Sterile-filter** the solution using a **0.22 µm vacuum filter unit**.
7. Aliquot if desired under sterile conditions.

### Labelling

Clearly label each bottle or aliquot with: - Buffer name: **ACK Lysis Buffer** - Date of preparation - Preparer initials - **Sterile - pH adjusted to 7.2**

### Storage and stability

- Store at **4 °C**.
- Stable for up to **3 months** when sterile.
- **Do not freeze**, as freezing may alter buffer performance.

### Safety and handling

- Wear gloves and eye protection; ammonium chloride is an irritant.
- Dispose of unused buffer and erythrocyte-containing waste according to biological/chemical waste regulations.
- For use in primary tissue processing, preparation and sterile filtration should be performed using appropriate aseptic technique.

## Version history

Version	Date	Author	Changes
v1.1	2026-01-09	Dillon Corvino	Corrected concentrations, explicit pH adjustment, clarified sterility and labelling requirements
v1.0	2025-11-21	Dillon Corvino	Initial Quarto buffer document for ACK lysis buffer