

33% Percoll Solution

Abdullah Lab, IMMEI, University Hospital Bonn

2025-11-20

Protocol ID: BUF-PERC-33-v1.0

Version: v1.0

Author: Dillon Corvino

Purpose

This buffer is used as a 33% Percoll solution for density-based enrichment of leukocytes from digested liver. It separates hepatocytes (floating) from leukocyte- and erythrocyte-enriched pellets during centrifugation.

Table of contents

| | |
|--|----------|
| Purpose | 1 |
| Linked protocols | 1 |
| Composition (per 10 mL 33% Percoll) | 2 |
| Preparation | 2 |
| Storage and stability | 2 |
| Reagent details | 2 |
| Safety (brief) | 3 |
| Version history | 3 |

Linked protocols

This buffer is used in:

- **Mouse Liver Dissociation to Single-Cell Suspension – LIV-001 (v1.0)**

Composition (per 10 mL 33% Percoll)

Based on the original SOP:

| Component | Stock | Volume for 10 mL | Approx. final composition | Notes |
|-------------------|-------|------------------|---------------------------|-----------------|
| Percoll (sterile) | 100% | 3.33 mL | 33% (v/v) | Undiluted stock |
| PBS 1× | — | 6.29 mL | — | Isotonic medium |
| PBS 10× | 10× | 370.37 μ L | Adjusts osmolarity | Added last |

Preparation

1. Equilibrate **Percoll stock** to **room temperature** before preparation.
2. For each **10 mL** of 33% Percoll:
 1. In a sterile 15 mL or 50 mL tube, add **3.33 mL** Percoll (sterile).
 2. Add **6.29 mL** PBS 1×
 3. Add **370.37 μ L** PBS 10×
 4. Mix gently by inversion until homogeneous.
3. If preparing multiple aliquots (e.g. **10 mL per liver**), scale volumes accordingly.
4. Aliquot **10 mL** per labeled 15 mL tube for each liver sample.
5. Label tubes with:
 - “33% Percoll”
 - Date of preparation
 - Version (BUF-PERC-33-v1.0)

i Note

Maintain sterility by preparing the solution in a biosafety cabinet or sterile environment if it will contact viable cells used for culture or sensitive downstream applications.

Storage and stability

- **Working 33% Percoll solution** is generally prepared **fresh on the day of use**.
- Short-term storage (same day) at **room temperature** is acceptable.
- Do not freeze 33% Percoll; discard any unused working solution at the end of the day.
- Stock Percoll (100%) should be stored according to manufacturer’s instructions (typically at room temperature), protected from contamination.

Reagent details

| Component | Supplier | Cat# | Notes |
|----------------|----------|-------|---|
| Percoll (100%) | [TBD] | [TBD] | Sterile Percoll; density gradient medium |
| PBS 1× | [TBD] | [TBD] | Sterile, Ca ²⁺ /Mg ²⁺ -free recommended |
| PBS 10× | [TBD] | [TBD] | Sterile concentrate; used to adjust osmolarity |

Safety (brief)

- Percoll solutions should be handled with gloves and standard lab PPE.
- Dispose of Percoll-containing waste according to local chemical and biohazard regulations when mixed with biological material.
- Avoid ingestion and contact with eyes; consult the manufacturer's safety data sheet (SDS) for details.

Version history

| Version | Date | Author | Change summary |
|---------|------------|----------------|---------------------------|
| v1.0 | 2025-11-20 | Dillon Corvino | Initial buffer definition |