

# R10 Medium (RPMI + 10% FCS + Pen/Strep)

Abdullah Lab, IMMEI, University Hospital Bonn

**Buffer ID:** BUF-R10-MED-001

**Version:** v1.0

**Author:** Dillon Corvino

**Date:** 2025-11-21

## Purpose

R10 is complete RPMI 1640 medium supplemented with 10% heat-inactivated FCS and Penicillin/Streptomycin. It is used for PBMC washing after Ficoll separation, for quenching CellTrace dyes, and as a general-purpose lymphocyte culture medium.

## Working buffer composition

Component	Final concentration	Notes
RPMI 1640	—	Base medium
Fetal calf serum	10% (v/v)	Heat-inactivated
Penicillin/Streptomycin	1% (v/v)	100× stock → 1× final

## Preparation

### Stock solutions required

Reagent	Stock concentration	Notes
Fetal calf serum (FCS)	—	Heat-inactivate if required
Penicillin/Streptomycin	100×	Add to 1× final
RPMI 1640	—	Sterile-filtered

## Preparation of 500 mL R10

1. Aliquot **450 mL** sterile RPMI 1640 into a sterile bottle.
2. Add **50 mL** FCS (10% final).
3. Add **5 mL** Penicillin/Streptomycin (100× stock) → **1% final**.
4. Gently invert to mix; **do not vortex**.
5. Label with:
  - *R10 complete medium*
  - Date of preparation
  - Expiry (4 weeks after preparation)



### Tip

Pre-warming R10 to 37 °C before use improves cell recovery following Ficoll separation or staining.

## Storage and stability

- Store at **4 °C**, protected from light.
- Use within **4 weeks** after preparation.
- Do not freeze.

## Reagent details

Reagent	Supplier	Cat. #	Notes
RPMI 1640	Various	—	Base medium
Fetal calf serum	Various	—	Heat-inactivated
Pen/Strep (100×)	Various	—	Final 1× concentration

## Safety

- Handle medium aseptically in a biosafety cabinet.
- Dispose of contaminated or expired medium as biohazardous waste.
- Follow institutional guidelines for handling FCS-containing media.

## Version history

Version	Date	Author	Changes
v1.0	2025-11-21	Dillon Corvino	Initial R10 buffer document.