

# LCMV Plaque Assay on Vero Cells

Abdullah Lab, IMMEI, University Hospital Bonn

2025-12-01

**Protocol ID:** MUS-LCMV-PLA-001

**Version:** v1.0

**Author:** Dillon Corvino

## Purpose

This protocol describes a standard plaque assay for determining the infectious titer of LCMV stocks using Vero cell monolayers. The readout is plaque-forming units per milliliter (PFU/mL).

### LCMV biosafety and handling

LCMV is a zoonotic arenavirus capable of causing severe disease in humans, including meningitis and congenital infection. Work with live LCMV must:

- Be conducted under the appropriate institutional **biosafety level (typically BSL-2 with additional local conditions)** and approved animal and biosafety protocols.
- Use **dedicated infectious CO incubators** and clearly labelled containment for all cultures and waste.
- Include suitable **PPE** (lab coat/gown, gloves, eye protection as appropriate) and measures to avoid sharps injuries, aerosol formation, and spills.
- Exclude or take special precautions for **pregnant or immunocompromised personnel** according to institutional policies.

All waste (liquid and solid) must be inactivated according to local regulations (e.g. chemical disinfection and/or autoclaving) **before** disposal.

## Table of contents

|  |          |
|--|----------|
| <b>Purpose</b>                               | <b>1</b> |
| <b>Principle</b>                             | <b>3</b> |
| <b>Materials</b>                             | <b>3</b> |
| Cells and Media . . . . .                    | 3        |
| Virus . . . . .                              | 3        |
| Reagents and Consumables . . . . .           | 3        |
| Equipment . . . . .                          | 3        |
| <b>Preparation</b>                           | <b>4</b> |
| Vero Cell Seeding . . . . .                  | 4        |
| Overlay Medium . . . . .                     | 4        |
| <b>Virus Dilution Series</b>                 | <b>4</b> |
| <b>Infection of Vero Cell Monolayers</b>     | <b>5</b> |
| <b>Overlay Application</b>                   | <b>5</b> |
| <b>Fixation and Staining</b>                 | <b>5</b> |
| <b>Plaque Counting and Titer Calculation</b> | <b>5</b> |
| <b>Controls and Quality Checks</b>           | <b>6</b> |
| <b>Troubleshooting</b>                       | <b>6</b> |
| <b>Safety considerations</b>                 | <b>6</b> |
| <b>Version History</b>                       | <b>7</b> |

## Principle

- Serial 10-fold dilutions of virus are prepared in serum-containing medium.
- Vero cell monolayers are infected with each dilution for a defined adsorption period.
- Inoculum is removed and cells are overlaid with a semi-solid medium to restrict virus spread to neighboring cells.
- After several days of incubation, individual plaques (zones of lysis/infection) are fixed, stained, and counted.
- Viral titer is calculated from plaque counts, dilution factor, and inoculum volume.

## Materials

### Cells and Media

- Vero cells (African green monkey kidney)
- Vero cell propagation medium (see `buf_vero_propagation_medium.qmd`)
- LCMV plaque assay medium (2× MEM-based overlay; see `buf_lcmv_plaque_assay_medium.qmd`)
- Methylcellulose or agarose for overlay (if preparing a gel-based overlay)

### Virus

- LCMV stocks to be titrated (Armstrong, Clone 13, WE, Docile, etc.)
- Known reference virus stock (optional internal control)

### Reagents and Consumables

- PBS without  $\text{Ca}^{2+}$  /  $\text{Mg}^{2+}$
- Crystal violet staining solution (e.g. 0.5% crystal violet in 20% ethanol)
- Fixative:
  - 4% paraformaldehyde in PBS, or
  - 10% neutral buffered formalin
- Sterile 6-well plates or 12-well plates
- Sterile conical tubes (15 mL, 50 mL)
- Micropipettes and sterile filtered tips
- Serological pipettes
- Waste containers with disinfectant (e.g. 10% bleach)

### Equipment

- Class II biosafety cabinet
- $\text{CO}_2$  incubator set to 37°C, 5%  $\text{CO}_2$
- Inverted light microscope
- Water bath (for overlay preparation if using agarose)
- Aspirator or vacuum system with appropriate trap and disinfectant
- Laboratory timer

## Preparation

### Vero Cell Seeding

1. One day before infection, seed Vero cells in 6-well plates to achieve a **90–100% confluent monolayer** on the day of infection.
2. Use Vero propagation medium and incubate at 37°C, 5% CO<sub>2</sub>.

#### Tip

Aim for evenly distributed, healthy monolayers. Avoid overconfluency and large areas of cell-free plastic.

### Overlay Medium

You can use either:

- **Methylcellulose-based overlay**, or
- **Low-melting agarose overlay**

A common approach is to mix **equal volumes** of:

- 2× plaque assay medium (see buffer file) and
- 2% methylcellulose or 2% agarose (cooled to ~42–45°C for agarose)

to obtain a final 1× medium with 1% methylcellulose or 1% agarose.

Prepare and keep overlay medium **warm (37–42°C)** until use.

### Virus Dilution Series

1. Thaw a virus aliquot rapidly in a 37°C water bath and immediately place on ice.
2. Prepare a series of 10-fold dilutions in Vero propagation medium (e.g. 10<sup>-1</sup> to 10<sup>-6</sup>).
3. Use sterile 1.5 mL or 15 mL tubes and mix gently by pipetting.

Example:

- Tube 1: 900 µL medium + 100 µL virus stock (10<sup>-1</sup>)
- Tube 2: 900 µL medium + 100 µL from Tube 1 (10<sup>-2</sup>)
- Continue to the desired dilution (e.g. 10<sup>-6</sup>).

## Infection of Vero Cell Monolayers

1. Remove growth medium from each well of the 6-well plate.
2. Rinse wells once with PBS (optional but recommended).
3. Add virus dilutions:
  - Typically 200–300  $\mu\text{L}$  per well.
  - Use at least duplicate wells per dilution if possible.
4. Gently rock plates to distribute inoculum evenly.
5. Incubate at 37°C, 5% CO<sub>2</sub> for **1 hour**.
6. Gently rock plates every 10–15 minutes during adsorption.

## Overlay Application

1. After the adsorption period, carefully aspirate the virus inoculum from each well without disturbing the cell monolayer.
2. Add **2–3 mL** of warm overlay medium to each well.
3. Swirl plates gently to ensure even coverage.
4. Incubate plates at 37°C, 5% CO<sub>2</sub> for **3–5 days**, depending on strain and desired plaque size.
  - Armstrong and WE often show clear plaques by day 3–4.
  - Clone 13 and Docile may require up to day 5 for optimal plaque definition.

## Fixation and Staining

1. At the chosen endpoint, remove plates from the incubator.
2. Carefully add an equal volume of **fixative** (e.g. 4% paraformaldehyde) directly on top of the overlay to fix cells and inactivate virus.
  - Incubate at room temperature for **30–60 minutes**.
3. Aspirate overlay and fixative into a disinfectant-containing waste container (e.g. 10% bleach).
4. Gently rinse wells with tap water or PBS.
5. Add sufficient **crystal violet staining solution** to cover the monolayer (e.g. 1–2 mL per well).
6. Incubate at room temperature for **15–30 minutes**.
7. Rinse thoroughly with tap water until excess stain is removed.
8. Air-dry plates.

Plaques will appear as clear (unstained) zones against a darkly stained monolayer.

## Plaque Counting and Titer Calculation

1. Count plaques in wells with **10–100 plaques** for accurate quantification.
2. For each counted well, record:
  - Dilution factor (e.g. 10<sup>-5</sup>)
  - Number of plaques
  - Inoculum volume (e.g. 0.2 mL)

### 3. Calculate PFU/mL:

$$\text{PFU/mL} = (\text{Number of plaques} \times \text{Dilution factor}) / \text{Volume of inoculum (mL)}$$

Example:

- 45 plaques at a  $10^{-5}$  dilution with 0.2 mL inoculum

PFU/mL is calculated as:

$$\text{PFU/mL} = (45 \times 10^5) / 0.2 = 2.25 \times 10^7 \text{ PFU/mL}$$

### 4. Average PFU/mL from replicate wells at the same dilution.

## Controls and Quality Checks

- Include at least one **mock-infected well** (no virus) to check for contamination and monolayer health.
- Optionally include a **reference virus stock** of known titer to verify assay consistency.
- Plaques should be well separated and have clearly defined borders.
- Monolayers should remain intact in mock-infected wells, with no evidence of cytotoxicity from medium or overlay.

## Troubleshooting

| Issue                        | Possible cause                                | Suggested action   |
|------------------------------|---|--|
| No plaques at any dilution   | Virus stock inactive or too dilute            | Check storage history; test a less diluted stock         |
| Too many plaques (confluent) | Dilutions not high enough                     | Extend dilution series (e.g. to $10^{-6}$ or $10^{-7}$ ) |
| Irregular or fuzzy plaques   | Overlay too thin or too thick; incubation off | Adjust overlay volume; verify incubator conditions       |
| Poor cell monolayer          | Cells not healthy or not fully attached       | Improve cell culture conditions before assay             |
| High background staining     | Incomplete washing after staining             | Increase number and volume of wash steps                 |

## Safety considerations

- Confirm that all work with LCMV and infected animals/samples is covered by approved **biosafety and animal use protocols**.
- Perform all live-virus manipulations (thawing virus stocks, infection steps, handling of unfixed plates) in a **certified Class II biosafety cabinet**.
- Use **dedicated infectious incubators** and clearly label plates, flasks, and waste containers with the virus strain and date.
- Wear appropriate **PPE** at all times (lab coat/gown, gloves, eye/face protection as required).

- Decontaminate all liquid waste (e.g. with freshly prepared bleach) before disposal. Solid waste (tips, plates, tubes) must be autoclaved or otherwise inactivated according to institutional guidelines.
- OPD and its solutions are handled as **chemical hazardous waste**; follow chemical safety rules in addition to biosafety rules.
- Pregnant or immunocompromised personnel should not work with LCMV unless explicitly allowed under institutional policy with additional safeguards.

## Version History

| Version | Date       | Description                    |
|---------|------------|--------------------------------|
| v1.0    | 2025-12-01 | Initial LCMV plaque assay SOP. |