

# NK Cell Culture Medium (cRPMI for Murine NK Cells)

Abdullah Lab, IMMEI, University Hospital Bonn

2025-11-21

**Protocol ID:** BUF-MUS-NK-MED-001

**Version:** v1.0

**Author:** Dillon Corvino

## Purpose

This medium is used for in vitro culture and proliferation assays of murine NK cells. It is based on complete RPMI 1640 supplemented with nutrients essential for lymphocyte viability and proliferation. **IL-15 is added fresh immediately before plating cells** and is not included in the stock medium.

## Working buffer composition

| Component                        | Final concentration | Notes                          |
|----------------------------------|---------------------|--------------------------------|
| RPMI 1640                        | —                   | Base medium                    |
| Fetal calf serum (FCS)           | 10% (v/v)           | Heat-inactivated               |
| Penicillin/Streptomycin          | 1% (v/v)            | 100× stock → 1× final          |
| HEPES                            | 10 mM ( 1% v/v)     | Buffering at 37 °C             |
| Sodium pyruvate (NaPy)           | 1 mM (1% v/v)       | Energy substrate               |
| Non-essential amino acids (NEAA) | 1% (v/v)            | Supports lymphocyte metabolism |
| -mercaptoethanol ( -ME)          | 50 M (e.g. 0.05 mM) | Reduces oxidative stress       |

### Note

IL-15 is **not included** in this stock medium. It must be added to a final concentration of 100 IU/mL immediately before culturing NK cells.

## Preparation

### Stock reagents required

| Reagent                | Stock concentration | Notes                 |
|------------------------|---------------------|-----------------------|
| RPMI 1640              | —                   | Sterile-filtered      |
| FCS, heat-inactivated  | —                   | 10% final             |
| Pen/Strep              | 100×                | 1% final              |
| HEPES buffer           | 1 M                 | 10 mM final ( 1% v/v) |
| Sodium pyruvate (NaPy) | 100 mM              | 1 mM final (1% v/v)   |
| NEAA                   | 100×                | 1% final              |
| -mercaptoethanol       | 50 mM or 14.3 M     | See dilution below    |

### Preparation of 500 mL NK medium (cRPMI base)

1. In a sterile bottle, combine:
  - **450 mL** RPMI 1640
  - **50 mL** FCS (10% final)
2. Add supplements:
  - **5 mL** Pen/Strep (100× stock → 1% final)
  - **5 mL** HEPES (1 M stock → 10 mM final)
  - **5 mL** Sodium pyruvate (100 mM stock → 1 mM final)
  - **5 mL** NEAA (100× stock → 1% final)
3. Add -mercaptoethanol:
  - Final concentration: **50 M**
  - If using 50 mM stock: add **0.5 mL**
  - If using 14.3 M stock: prepare a working dilution before adding
4. Mix gently by inversion; do **not vortex**.
5. Label with:
  - “NK medium (cRPMI base)”
  - Date of preparation
  - Expiry (4–5 weeks at 4 °C)

### Warning

-mercaptoethanol is hazardous. Handle it in a fume hood and avoid skin contact. Do not open concentrated -ME bottles inside a biosafety cabinet.

## Usage

- Warm medium to **37 °C** before adding cells.
- **Add IL-15 fresh** immediately prior to plating NK cells:
  - Final concentration **100 IU/mL**.
- Do not store IL-15 in the base medium.

## Storage and stability

- Store at **4 °C**, protected from light.
- Use within **4–5 weeks** of preparation.
- Do **not freeze**.

## Reagent details

| Reagent             | Supplier      | Cat. # | Notes                       |
|---------------------|---------------|--------|-----------------------------|
| RPMI 1640           | Various       | —      | Base medium                 |
| FCS                 | Various       | —      | Heat-inactivated            |
| Pen/Strep           | Various       | —      | 100× working stock          |
| HEPES               | Various       | —      | 1 M stock                   |
| Sodium pyruvate     | Various       | —      | 100 mM stock                |
| NEAA                | Various       | —      | 100× stock                  |
| -mercaptoethanol    | Sigma/etc     | —      | Handle with caution         |
| IL-15 (added fresh) | PeproTech/etc | —      | Not included in base medium |

## Safety

- Handle -mercaptoethanol and all supplements using gloves and eye protection.
- -ME requires fume hood handling due to volatility.

- Prepare medium under aseptic conditions inside a biosafety cabinet.
- Dispose of expired or contaminated medium as biohazardous waste.

## Version history

| Version | Date       | Author         | Changes                                      |
|---------|------------|----------------|--|
| v1.0    | 2025-11-21 | Dillon Corvino | Initial Quarto buffer document for NK medium |