

Standard Flow Cytometry Staining (Surface, Intracellular, Intranuclear)

Abdullah Lab, IMMEI, University Hospital Bonn

2025-12-29

Protocol ID: GEN-FLOW-STD-001

Version: v1.1

Author: Dillon Corvino

Purpose

This protocol describes a standard flow cytometry staining workflow for suspension cells, including:

- Optional **tetramer staining**.
- Surface staining with Fc block and live/dead dye.
- Optional intracellular staining (e.g. BD Cytfix/Cytoperm-based).
- Optional intranuclear staining (e.g. eBioscience FOXP3 fix/perm-based).

It is intended as a general framework for murine or human samples and can be adapted to specific panels and antibody combinations.

Buffer nomenclature (important)

- **FACS Buffer** = PBS + FCS
- **MACS Buffer** = PBS + FCS + EDTA

Unless otherwise specified, **MACS buffer is used throughout this protocol** to minimise cell clumping during staining and washing.

Tetramer staining is a specific exception (see below).

Critical notes (read before starting)

- Keep cells and antibody mixes cold (on ice or at 4 °C) and protected from light whenever possible.
- Default staining range: **$2 \times 10^5 - 1 \times 10^6$ cells per well** (96-well format).
- Always include appropriate controls (unstained, live/dead only, FMO, single-colour controls).
- Fixation and permeabilization conditions must be compatible with antibodies and fluorochromes used.
- Commercial fixation/permeabilization kits should be used according to manufacturer instructions.

Warning

Improper fixation/permeabilization or incompatible buffers can destroy epitopes and fluorochromes, leading to misleading intracellular or intranuclear staining.

Warning

Antibodies and tetramers are **protein reagents and mechanically sensitive**.

Do **not vortex**, as this can denature proteins and promote aggregation.

Instead, mix gently by pipetting or inversion. If needed, briefly spin reagents in a **mini “pulse” centrifuge (quick spin)** to collect liquid and pellet aggregates before use.

Tip

Surface markers (including live/dead dyes and tetramers) should generally be stained **before fixation/permeabilization** for best performance.

Approximate timing

- Tetramer staining (optional): 20 min
- Plate setup and Fc block: 10–15 min
- Surface staining: 30 min
- Fixation only (no intracellular): 15–20 min
- Fix + permeabilization + intracellular/intranuclear staining: 1.5–2 h
- Optional overnight intranuclear stain: +1 day

Total (surface only): ~1–1.5 h

Total (full workflow): ~3–4 h (or overnight)

Table of contents

Purpose	1
Buffer nomenclature (important)	1
Critical notes (read before starting)	2
Approximate timing	2
Procedure	4
Step 1 – Plate setup and buffers	4
Step 2 – Tetramer staining (optional, if applicable)	4
Step 3 – Surface staining (with live/dead)	4
3.1 If tetramer staining was performed (from Step 2)	5
3.2 If no tetramer staining was performed	5
3.3 Wash steps (both cases)	5
3.1 Surface staining only	5
Step 4 – Fixation and permeabilization	5
4.1 Initial fixation/permeabilization	6
4.2 Optional pause point	6
Step 5 – Intracellular or intranuclear staining	6
Materials	7
Reagents	7
Buffers used	7
Version history	7

Procedure

Step 1 – Plate setup and buffers

1. Prepare buffers:
 - **FACS buffer** (PBS + 2% FCS) — *for tetramer staining only.*
 - **MACS buffer** (PBS + 2% FCS + 2 mM EDTA) — *default buffer for all other steps.*
 2. Prepare single-cell suspensions.
 3. Count cells and adjust concentration to plate **$2 \times 10^5 - 1 \times 10^6$ cells per well.**
 4. Aliquot cells into a 96-well U- or V-bottom plate.
 5. Keep plates on ice or at 4 °C.
-

Step 2 – Tetramer staining (optional, if applicable)

This step applies when staining with MR1 or CD1d tetramers.

1. Prepare a **tetramer staining mix in FACS buffer (no EDTA)**:
 - Fc block: **1:100**
 - Tetramer(s): per manufacturer recommendation
 - Final volume: **50 µL per well**
2. Add **50 µL** tetramer mix to each well.
3. Incubate for **20 min on ice**, protected from light.

i Note

Tetramer staining is performed in **FACS buffer without EDTA** to avoid potential interference with multimer-TCR interactions.
Whether EDTA exclusion is strictly necessary remains under investigation and will be tested within the lab.

i Note

MR1-5-OP-RU tetramer can be added together with surface antibody staining.
For CD1d tetramers, compatibility with simultaneous surface staining is currently unclear and will be empirically tested.

4. Without washing, proceed directly to Step 3.
-

Step 3 – Surface staining (with live/dead)

Antibody titration guidance (applies to all surface staining): - Surface antibodies: default **1:200** if titration data are unavailable. - Transcription factors or weakly expressed markers: default **1:50**.

3.1 If tetramer staining was performed (from Step 2)

1. Prepare a **2× surface staining master mix in MACS buffer** containing:
 - Live/dead dye (2× final concentration)
 - Surface antibodies (2× final concentration)
2. Add **50 µL** of the 2× surface mix to each well already containing **50 µL tetramer mix**
→ **Final volume: 100 µL per well (1× final concentration)**.
3. Incubate for **30 min at 4 °C**, protected from light.

3.2 If no tetramer staining was performed

1. Prepare a **1× surface staining master mix in MACS buffer** containing:
 - Live/dead dye (1× final concentration)
 - Surface antibodies (1× final concentration)
2. Add **50 µL** of the 1× surface mix directly to each well.
3. Incubate for **30 min at 4 °C**, protected from light.

3.3 Wash steps (both cases)

1. Top up each well to **200 µL** with MACS buffer.
2. Centrifuge at **400 g for 2 min**.
3. Flick off the supernatant.
4. Wash once more with **200 µL MACS buffer**, centrifuge, and flick.

i Note

All antibodies should be titrated for optimal performance. Default dilutions are provided as a starting point only.

3.1 Surface staining only

1. Resuspend pellets in **100–200 µL 4% PFA in PBS**.
2. Store at **4 °C**, protected from light, until acquisition.

Step 4 – Fixation and permeabilization

This step assumes use of: - **BD Cytofix/Cytoperm** (cytoplasmic targets), or
- **eBioscience FOXP3 fix/perm** (intranuclear targets).

4.1 Initial fixation/permeabilization

1. Add **100 µL** fix/perm solution per well.
2. Incubate **20 min on ice**, protected from light.
3. Top up to **200 µL** with 1× PermWash.
4. Centrifuge and flick.
5. Resuspend in **200 µL** PermWash.
6. Centrifuge and flick.

4.2 Optional pause point

1. Resuspend pellets in **200 µL MACS buffer**.
2. Store at **4 °C overnight**, protected from light.

Note

This pause point is compatible with most intracellular workflows, but antigen stability should be confirmed for each marker.

Step 5 – Intracellular or intranuclear staining

1. Prepare an **intracellular or intranuclear antibody master mix** in the appropriate **PermWash** buffer (kit-specific):
 - Final volume: **50 µL per well**
2. Resuspend each cell pellet in **50 µL** of the antibody master mix.
3. Incubate **30 min on ice**, protected from light, **or** incubate **overnight at 4 °C**.

Tip

For **intranuclear targets** (e.g. **FOXP3, Helios, T-bet, Eomes**), overnight antibody incubation at **4 °C** can improve signal intensity, particularly for low-abundance or epitope-masked transcription factors.

Ensure plates are **well sealed, clearly labelled, and protected from light** during overnight incubation.

Warning

Overnight antibody incubation is **not recommended for cytoplasmic intracellular targets** (e.g. cytokines, granzymes, perforin). Prolonged incubation can increase background staining and does not typically improve signal for these antigens. Use **30–60 min** incubations for intracellular (non-nuclear) targets.

4. Top up each well to **200 µL** with PermWash.
5. Centrifuge and flick off the supernatant.
6. Wash once more with **200 µL PermWash**, centrifuge, and flick.

7. Resuspend each pellet in **100–200 μ L 4% PFA in PBS.**
8. Store at **4 °C**, protected from light, until acquisition.

Materials

Reagents

Reagent	Notes
PBS 1×	Ca ² /Mg ² -free
FCS (heat-inactivated)	2% v/v
EDTA	2 mM (MACS buffer only)
Fc block	Mouse: anti-CD16/32
Live/dead dye	As per manufacturer
Tetramers (e.g. MR1–5-OP-RU, CD1d)	Optional
Surface antibodies	Titrate individually
Intracellular/intranuclear antibodies	Kit-compatible
Fix/perm kits	BD or eBioscience
4% PFA	Fixation

Buffers used

- **FACS buffer** – PBS + 2% FCS (tetramers only)
- **MACS buffer** – PBS + 2% FCS + 2 mM EDTA (default)

Version history

Version	Date	Author	Change summary
v1.0	2025-11-21	Dillon Corvino	Initial version
v1.1	2025-12-29	Dillon Corvino	Added tetramer workflow, clarified buffers, titration defaults, overnight intranuclear staining